



# AP *medical*

## CS *line*



Organizzazione con Sistema di Gestione certificato  
Company with Management System certified

ISO 9001:2000

**SINERT**



Disposable and unpiercable square containers with turning lid for disposable of sharp and cutting waste

# CS line

Wide opening with indents for needle's removal

Temporary closure: one release  
Definitive safety lock: two releases



**CS 2 LT**

**CS 3 LT**

**CS 4 LT**

**CS 5 LT B**

**CS 5 LT H**

**CS 6 LT H**



- Easy and practical closure
- Wide range with 5 different sizes
- Wide opening to easy the discharging operation
- Device for insulin needles removal
- 5lt and 6 lt suitable for long needles disposal
- During incineration no toxic or harmful gas emission occurs
- Filling level indicated
- Table, wall and trolley supports available
- Autoclavable at 134° C for 25 minutes



BS 7320:1990

Item	Lt	Max load	Product cm	Weight kg	Pcs/carton	Carton cm	Pcs/pall	Pallet cm
CS 2 LT	2	1,30	14,5X14,5X19,3	0,24	60	79X40X34	1080	80X120X220
CS 3 LT	3	1,70	16,7X16,7X22,5	0,30	50	79X40X34	900	80X120X220
CS 4 LT	4	2,00	16,7X16,7X28,5	0,35	40	79X40X34	720	80X120X220
CS 5 LT B	5	2,60	20,2X20,2X20,0	0,36	30	79X40X34	540	80X120X220
CS 5 LT H	5	2,80	16,7X16,7X37,2	0,41	40	79X40X42	600	80X120X230
CS 6 LT H	6	3,40	16,7X16,7X48,2	0,46	30	79X40X34	540	80X120X220

# Instructions of use

## 1 Lid assembly




- Press firmly on the opposite corners



- Press firmly on the opposite sides

## 2 Needle removal



- Respect the maximum filling level indicated on the container

## 3 Lid closure



- **1 click**  
Temporary closure

- **2 click**  
Definitive safety lock


# Supports

- Table supports with suction cups available for entire CS line



- Supports for vertical trolley fixing available for the entire PBS line



 **AP** medical

Retailer's stamp

